




Department of
Primary Industries and
Regional Development

Protect
Grow
Innovate

An introduction to Nanopore DNA sequencing and data analysis

Asad Prodhan

Aims

- ❑ DNA extraction to Data generation
 - ❑ Some key points in each step
 - ❑ Now I have my data, how do I get started?
 - ❑ Some resources
- 

DNA extraction to Nanopore sequencing: an example



MinION
High Throughput Sequencing



Bacterial IDs

DNA extraction to Nanopore sequencing: an example



DNA
Extraction



Library
Preparation



MinION

High Throughput Sequencing



Sequencing
Reads

Read 1: ATGCCGGTATTAA ...

Read 2: CGGATTAACGAGTAT ...

hundreds of thousands of reads

100-500 GB per run




Bioinformatics
analysis




Bacterial IDs

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



DNA Extraction




https://lab.loman.net/2018/05/25/dna-extraction-book-chapter/#__RefHeading__Toc505877552

- High Molecular Weight (HMW) DNA = ~ 20-100 Kb
- longer fragments = longer reads
- longer reads are suitable for *De novo* assembly and structural variations (deletion, insertion, duplication, inversion, SNP, translocation)
- Pipet slowly to minimise shearing DNA
- QC is required
- DNA fragment size selection is possible.

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



DNA Extraction → Library Preparation



MinION
High Throughput Sequencing

Sequencing Reads

Read 1: ATGCCGGTATTAA ...
Read 2: CGGATTAACGAGTAT ...
*hundreds of thousands of reads
100-500 GB per run*

Bioinformatics analysis

Bacterial IDs

dsDNA: μg to pmol

DNA Length (bp)

1000

DNA Amount (μg)

0.066


0.1pmols of DNA

Clear

Calculate


- Take the recommended amount of input DNA
- How do we calculate fmol?
 - <https://www.promega.com.au/resources/tools/biomath/>
 - 1 pmol = 1000 fmol
 - 100 fmol = 1 Kb X 66 ng
- Prepare the ends of the DNA fragments
- Attach the barcodes and adapter sequences (need to be trimmed)

Library Preparation



Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



MinION




Fig. MinION Mk1C

- Run time parameters
- Basecalling mode: fast, hac (high accuracy) or sup (supper accuracy)
- Min qscore = quality threshold for the reads to pass
- Minimum length filter = length threshold for the reads to pass
- Data output format: fast5 or pod5 (before basecalling), fastq (after basecalling)

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



⇒ DNA Extraction ⇒ Library Preparation



MinION
High Throughput Sequencing



Sequencing Reads
Read 1: ATGCCGGTATTAA ...
Read 2: CGGATTAACGAGTAT ...
*hundreds of thousands of reads
100-500 GB per run*



Bioinformatics analysis

⇒ Bacterial IDs

Sequencing Reads


- fast5 & fastq format
- SFTP



- Transfer data from MinION to your Linux computer
- `sftp minit@MinION_Number.local:/data; then get -r *`
- fast5 will allow you to basecall again later if needed
(<https://github.com/asadprodhan/GPU-accelerated-guppy-basecalling>)
- fast5 takes a lot of storage space
- storage space can be reduced by converting fast5 to slow5
(<https://github.com/asadprodhan/Reducing-storage-space-by-converting-Nanopore-fast5-to-slow5-using-slow5tools>)

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



Bioinformatics Analysis


- Read relevant papers
- Identify bioinformatics tools
- How to install the tools?



you get an answer

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



Bioinformatics Analysis

Step 1: Install conda

- ☐ How to check if conda already installed
- ☐ If not, how to install conda in your computer



Step 2: Create conda environment

- ☐ Why conda environment
- ☐ How to create, use, and share conda environment



Step 3: Install softwares within your conda environment

- ☐ How to install softwares using conda



Step 4: Perform your bioinformatics analysis

- ☐ How to download sequences from NCBI
- ☐ How to perform BLAST search
- ☐ How to sort out the BLAST hits
- ☐ How to extract sequences of the BLAST hits


<https://github.com/asadprodhan/A-beginner-s-guide-to-Bioinformatics>



you get an answer

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



DNA
Extraction



Library
Preparation



MinION
High Throughput Sequencing



Sequencing
Reads

Read 1: ATGCCGGTATTAA ...
Read 2: CGGATTAACGAGTAT ...
*hundreds of thousands of reads
100-500 GB per run*



Bioinformatics
analysis



Bacterial IDs

Bioinformatics Analysis





Fig. Bioinformatics analysis workflow



you get an answer

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



Bioinformatics Analysis





Fig. Bioinformatics analysis workflow

you get an answer

- Genome assembly
(<https://github.com/asadprodhan/Setting-up-and-running-Miniasm-assembler>)
- Blastn
(<https://github.com/asadprodhan/A-beginner-s-guide-to-Bioinformatics>)
- Kraken2
(<https://github.com/asadprodhan/How-to-run-Kraken2-on-HPC-using-Singularity-container-and-Nextflow>)
- ANI
(<https://github.com/asadprodhan/Average-Nucleotide-Identity-ANI-analysis>)
- UBCG & Phylogeny
(<https://github.com/asadprodhan/Bacterial-phylogenetic-tree-reconstruction-using-UBCG-pipeline>)

Dissecting Each Step

DNA extraction to Nanopore sequencing: an example



DNA Extraction → Library Preparation



MinION
High Throughput Sequencing

Sequencing Reads

Read 1: ATGCCGGTATTAA ...
Read 2: CGGATTAACGAGTAT ...
hundreds of thousands of reads
100-500 GB per run

Bioinformatics analysis

Bacterial IDs

Bioinformatics Analysis





Fig. Bioinformatics analysis workflow


HPC Cluster
Nextflow & Singularity container
A collection of tutorials to learn coding from basics to bioinformatics.

(<https://www.linkedin.com/in/asadprodhan/>)



Summary

Analysis Steps



Example tools

NanoPlot

NanoFilt

NanoPlot

Miniasm, Flye, Canu, Spades etc

Medaka, NanoPolish

Quast

Blastn, Kraken2, ANI, Phylogeny etc

Acknowledgements

- Monica Kehoe
- Tonny Kinene
- Dominie Wright
- Cuiping Wang



Department of
**Primary Industries and
Regional Development**



Thank you

dpird.wa.gov.au    

Important disclaimer

The Chief Executive Officer of the Department of Primary Industries and Regional Development and the State of Western Australia accept no liability whatsoever by reason of negligence or otherwise arising from the use or release of this information or any part of it.

Copyright © State of Western Australia (Department of Primary Industries and Regional Development), 2022.